TOTAL **HOSPITAL DISINFECTANT/DEODORANT**



DESCRIPTION

This disinfectant spray is ideal for all surfaces. It is more than just an air spray. Its phenylphenol base and country garden and fresh & clean scents make it the perfect product for total rest room sanitation. This product is effective against a wide variety of microorganisms including staph, salmonella, pseudomonas, HIV-1 (AIDS virus), TB (tubercle bacilli), Streptococcus pyogenes, Influenza A, (Hong Kong), and Herpes simplex I & II. It is also fungicidal against pathogenic fungi and prevents the growth of mold and mildew. Furthermore, its residual, selfsanitizing activity controls germs for up to 14 days on fabric and up to 21 days on hard surfaces after use! In addition, it eliminates airborne malodors from sick rooms and musty closets. Also, it does not contain any ozone-depleting CFC's.

* Schools

* Institutions

* Nursing Homes

* Shopping Malls

* Show Rooms

* Country Clubs

* Restaurants

* Cafeterias

* Hospitals

BENEFITS - FEATURES

- * Kills Influenza A₂ (Hong Kong)
- * Kills HIV-1 (AIDS)
- * Kills Herpes simplex I & II
- * Kills TB
- * Deodorizes * Disinfects
- * No CFC's

* Residual

ADVANTAGES

- * Self-sanitizing action controls germs on fabric for up to 14 days and on hard surfaces for up to 21 days
- * Eliminates all airborne malodors at their source including odor-causing bacteria
- RECOMMENDED FOR
 - * Office Buildings
 - * Government Buildings
 - * Hotels
 - Motels
 - Churches
 - * Airports
 - * Factories
 - * Supermarkets

DIRECTIONS

WHERE TO USE

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

For Precleaning: Preclean surfaces with a suitable detergent prior to disinfecting.

For Disinfecting: Hold dispenser upright 6 to 8 inches from surface. Spray 2 to 3 seconds or until entire surface is wet. Allow to remain at least 10 minutes. Allow to air-dry or remove with clean, dry cloth.

For Deodorizing: Hold dispenser upright. Spray toward center of room 3 to 5 seconds.

For Mildew Control: Spray surface until moist. Do not use on polished, wood furniture or rayon fabrics. Repeat application periodically to maintain control.

*KILLS HIV-1 ON PRECLEANED, ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/ BODY FLUIDS in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/ objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION

AGAINST HIV1 OF SURFACES/ **OBJECTS SOILED WITH BLOOD/BODY**

FLUIDS: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks, and/or eve coverings

CLEANING PROCEDURES: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application.

CONTACT TIME: Allow surface to remain wet for 10 minutes.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

SPECIFICATIONS

Spray pattern.....Ultra-fine mist Appearance.....Clear, light yellow Fragrance.....Fresh & CleanCountry Garden

Boiling point (Concentrate)	176°F
Solubility in water	Complete
Specific gravity (H ₂ 0=1)	0.83
рН	

Vapor density (Air=1)	55 to 65
Flammability	Flammable
Storage stability	130°F max
Net weight	16 oz. (454 g)

Another Exclusive Product of





١

* Liaht Swithes Bedpans

Bed Springs & Frames

Restroom Fixtures

- * Diaper Pails
- * Toilet Seats
- * Wheelchairs
- * Garbage Cans * Flush Handles * Telephones

* All Airborne Odors

* Door Knobs

* Linen Carts

* Hampers

* Walls

* Sinks